

Acid Pump/Reagent Dispenser

Catalog No. F11656-0000



The Reagent Dispenser comes with a length of polyethylene tubing 1/4" O.D. Insert the end of this tube into the 1/4" hole on the underside of the body of the dispenser. Cut the tube to the length needed for the size of bottle being used. The tube should reach the bottom of the bottle. If the tube becomes loose in the hole, expand the end by putting it into very hot water, let cool and reinsert.

A short length of 3/8" diameter polyethylene tube is supplied for the outlet. Insert this tube into the side hole of the dispenser body.

For use on a 5 pint reagent bottle, place the 1/8" thick tapered washer on the bottle top and screw on the dispenser until tight.

Do not overtighten.

For use on a 1/2 gallon or 1 gallon bottle, use the 1/2" wide sleeve in the same manner.

Screw the small bottle into the threaded hole in the side of the dispenser body.

To achieve maximum delivery with minimum effort use rapid, light finger tip squeezes on the small bottle.

Fluid will be dispensed from the short tube. Light pressure on the relief valve button on top of the dispenser stops flow quickly.

You May Also be Interested in:

Cool Stir® Large Volume Magnetic Stirrer

15½ x 15½ x 9½" (39 x 39 x 23cm)

The Cool Stir® Large Volume Magnetic Stirrer efficiently stirs fluid volumes from 10 liters up to 208 liters (2.5 to 55 gallons). Stirring speed is continuously changeable to around 600 rpm (maximum operational rpm depends on volume and density) delivering torque for mixing through powerful agitation with a deep vortex.

- Direct drive motor rotates a powerful rare earth magnet that couples optimally with 3" through 6" stirring bars (F37122 series and F37118 series- not included)
- Rugged built-in aluminum support system and durable polypropylene platform are engineered to carry up to 220kg (485lbs); permits mixing in a regular plastic barrel or drum
- Unit and its 100 watt (1/8 horsepower) motor are air cooled, and heat transmission to thermolabile solutions is prevented by the use of insulation
- Polypropylene sides are simple to clean
- DC motor: 100 watts, 50/60 Hz, 3-wire grounded cord
- Net weight: 15 kg (33lbs.); shipping weight: 17.3 kg (38lbs.); ships oversize

Catalog No. F37028-0001



SP Bel-Art

2024 Broad Street, Pocomoke, MD 21851 USA

Tel: 973-694-0500 | cservice@scientificproducts.com

belart.com

Individual use varies. SP Bel-Art assumes no obligation or liability for any advice furnished by it or for results obtained with respect to these products. All advice is given and accepted at the user's risk. For warranty information and to register your product visit scientificproducts.com.

