

Clavies® Biohazard Autoclave Gloves

Catalog No. H13201-0000, H13201-0001

**UK
CA**

Clavies® Gloves are suitable for temperatures up to 232°C (450°F). Always exercise caution when working with hot objects and equipment. Contact time with hot objects should always be minimized.

DO NOT USE A WET GLOVE

DO NOT HOLD GLOVE IN LIVE STEAM

Clavies® Gloves will absorb liquids that conduct heat energy, which may cause burns to the hand.

Be careful when picking up sharp or pointed objects that could penetrate the glove.

Clavies® Gloves are made of cotton terry cloth. Hand washing and air drying is recommended.



SP Bel-Art

2024 Broad Street, Pocomoke, MD 21851 USA

Tel: 973-694-0500 | cservice@scientificproducts.com

belart.com



Individual use varies. SP Bel-Art assumes no obligation or liability for any advice furnished by it or for results obtained with respect to these products. All advice is given and accepted at the user's risk. For warranty information and to register your product visit scientificproducts.com.

