

Poly-Jaque™ Plastic Lab Jack

Catalog No. F18391-0001

1. Using the control knob, check the operation of the jack for its full range up and down. Although the jack will support up to 15 pounds (6.8kg.), do not attempt to lift the weight from the closed position.
2. Raise the platform to a level slightly higher than needed and then place the weight centrally on the platform.
3. Minor adjustments can now be made with the control knob. The jack can be placed in a solution bath or corrosive atmosphere without fear of rusting, corroding or conducting electricity.

Closed height of the jack is 3" (76mm), maximum open height is 12" (300mm).

The jack can be placed in a steam autoclave in an unloaded (unstressed) condition at 250°F (121°C).

The bottom plate of the jack may be drilled for permanently attaching the jack to a support.



SP Bel-Art

2024 Broad Street, Pocomoke, MD 21851 USA

Tel: 973-694-0500 | cservice@scientificproducts.com

belart.com

Individual use varies. SP Bel-Art assumes no obligation or liability for any advice furnished by it or for results obtained with respect to these products. All advice is given and accepted at the user's risk. For warranty information and to register your product visit scientificproducts.com.

