

Lab-Aire® II Electric Labware Drying Rack

Catalog No. F18819-0013 • F18819-0023

Controls:

- Main Power Switch: Push-button with integral pilot light.
- Heat Control: Stepless time proportional controller.
- Circuit Breaker: 250VAC/50VDC, 8A

Operation:

Turn the main power switch and heat control to OFF. Plug line cord into an appropriate outlet. Turn main power switch to ON; pilot light should come on. Adjust heat control as desired. Invert glassware and place on pegs.

Pegs & Plugs:

On each panel, pegs should be arranged from small on top rows to large on bottom rows. On multipanel units, position small pegs on upper panels and large pegs on lower panels. Plugs are provided to cover holes where no pegs are desired.

Baskets are provided to hold spare unused pegs and plugs.

Maintenance:

Dust filter, which is mounted at the underside of the base, should be cleaned once a month (115 hours of operation) by vacuuming. The filter is removed by snapping off its holder.



Glassware Dryer Accessory Pegs - Multi-pack

Description	Catalog No.
includes: 19 small pegs, 19 medium pegs, 19 large pegs, and 6 plugs	C18931-0004



SP Bel-Art

2024 Broad Street, Pocomoke, MD 21851 USA

Tel: 973-694-0500 | cservice@scientificproducts.com

belart.com

Individual use varies. SP Bel-Art assumes no obligation or liability for any advice furnished by it or for results obtained with respect to these products. All advice is given and accepted at the user's risk. For warranty information and to register your product visit scientificproducts.com.

