

Safety Goggles

Catalog No. H24873-0000

Please Read This Before Wearing

WARNING: SP Bel-Art Safety Goggles offer limited protection against moderate impact hazards, metal chips, particles, and sparks. These lenses are not unbreakable or impenetrable. This polycarbonate product offers greater impact resistance than glass lenses, and meets the performance standards of ANSI Z87.1-2003. Please refer to your ANSI manual or safety specialist for proper selection as some applications may require face-shields, machine guards, etc. in addition to these goggles. If the lenses become scratched or pitted, this product should be discarded and replaced immediately.

Cleaning: Use the SP Bel-Art LabLens® Liquid (F24833-0000) or the SP Bel-Art Disposable Lens Cleaning Station (F24835-0000) to clean goggles. Mild soap and warm water may also be used. DO NOT use alcohol, solvents, or special cleaners. Air dry or use a soft tissue to pat dry the lenses.



SP Bel-Art

2024 Broad Street, Pocomoke, MD 21851 USA

Tel: 973-694-0500 | cservice@scientificproducts.com

belart.com

Individual use varies. SP Bel-Art assumes no obligation or liability for any advice furnished by it or for results obtained with respect to these products. All advice is given and accepted at the user's risk. For warranty information and to register your product visit scientificproducts.com.

