

Mercury Collector Replacement Pads

Catalog No. F36251-0000

These open cell polyurethane pads are replacements for the pad supplied with the Bel-Art mercury collector F32650-0000. When handling the used pad, wear gloves and a mask to prevent contact with hazardous mercury. Dispose of the used pad in accordance with local regulations.

1. Using a scraper, such as a putty knife, remove the used pad from the three legged platform and scrape off excess dried cement.
2. To install a new pad on the platform, use rubber cement (usually available at art supply stores). Apply a layer of cement to the platform and center the pad on the platform with light pressure. For your safety, comply with the warnings on the rubber cement instructions.

**SP Bel-Art**

2024 Broad Street, Pocomoke, MD 21851 USA

Tel: 973-694-0500 | cservice@scientificproducts.com

belart.com

Individual use varies. SP Bel-Art assumes no obligation or liability for any advice furnished by it or for results obtained with respect to these products. All advice is given and accepted at the user's risk. For warranty information and to register your product visit scientificproducts.com.



936251004 Rev. 2