

I N S T R U C T I O N S



Bel-Art Bulb Pipettor

Catalog No. H37888-0000

Hold pipette as usual and hold the Red Bulb Pipette in the other hand with the open end down. Squeeze the bulb and touch the open end over the top of the pipette. Release the squeeze on the bulb and allow the fluid to rise in the pipette. Quickly remove the bulb and stop fluid outflow with your index finger as usual.

***May be autoclaved after disassembly.
Allow parts to cool before reassembly.***

If any solution is aspirated inside bulb, take apart bulb and nose, wash, dry and reassemble by twisting and pushing nose and bulb together.



Help others choose the right product. Write a review at www.belart.com

661 Route 23 South, Wayne, NJ 07470

TEL: 973-694-0500 • FAX: 973-694-7199 • www.belart.com

Bel-Art – SP Scienceware assumes no obligation or liability for any advice furnished by it, or for results obtained with respect to these products.

All such advice is given and accepted at the buyer's risk.

© 2020. All Rights Reserved.

Bel-Art® – SP Scienceware® is a division of SP Industries, Inc.



937888001 - Rev. 2

