

Mini-Vacuum Desiccator

Catalog No. 41099-0000

The Mini-Vacuum Desiccator consists of a clear top section, a white bottom section, an O-ring and a vacuum retention valve.

The vacuum retention valve fits on the tapered neck at the top. It accepts 1/4" inside diameter tubing.

Place the O-ring in the half round groove in the bottom plate.

Insert the sample to be desiccated.

Place the clear desiccator top onto the base, making sure that the bottom edge comes into contact with

the O-ring. Press down until the flange snaps into the two locks.

To remove the top, turn the valve to shut off the vacuum and then turn to vacuum release position. Press back one of the snap locks and raise the top section.



You May Also be Interested in:

Secador Desiccator Cabinets

Demand More From Your Desiccator

Select the Secador™ Desiccator Cabinet that suits your lab!

- Patented sealed construction – Airtight, dust and moisture free storage
- Heavy duty copolyester composition builds rugged cabinet with long service life
- Imbedded UV protection blocks 99% of light from cabinet and contents
- Superior Chemical Resistance- Prevent staining and crazing
- Built in Hygrometer for easy RH monitoring
- Stackable – Saves valuable bench space
- Use with any Scienceware® desiccant



Secador™ 1.0
Volume .75 cubic feet



Secador™ 4.0 Horizontal & Vertical
Volume 2.0 cubic feet

Description	Catalog Number
Secador™ 1.0 Desiccator Cabinet	42071-0000
Secador™ 4.0 Desiccator Horizontal Cabinet	42074-0000
Secador™ 4.0 Desiccator Vertical Cabinet	42074-1000

Your opinion counts! Write a product review at www.sp-belart.com



Bel-Art Products, a part of SP

661 Route 23 South, Wayne, NJ 07470 USA

Tel: 973-694-0500 | cservice.belart@spindustries.com

sp-belart.com

Individual use varies. Bel-Art Products/SP assumes no obligation or liability for any advice furnished by it or for results obtained with respect to these products. All advice is given and accepted at the user's risk. For warranty information visit belart.com/corporate. For warranty information and to register your product visit sp-belart.com/corporate

