





Spinbar® Magnetic Sink Strainer

Catalog No. F37787-0000 (blue) • F37787-0001 (green)



INSTRUCTIONS:

- Place strainer in sink drain in an upright position so that Scienceware® logo is visible.
- Flat edge allows strainer to be placed close to sink walls.
- Four strong magnets in each corner attract loose Spinbar® magnetic stirring bars and prevent them from falling into the drain.
- Periodically inspect and clean strainer as needed.



Help others choose the right product. Write a review at www.belart.com

661 Route 23 South, Wayne, NJ 07470

TEL: 973-694-0500 • FAX: 973-694-7199 • www.belart.com

937902113 – Rev. 2


Bel-Art – SP Scienceware assumes no obligation or liability for any advice furnished by it, or for results obtained with respect to these products. All such advice is given and accepted at the buyer's risk. © 2019. All Rights Reserved.

Bel-Art® – SP Scienceware® is a division of SP Industries, Inc.




Spinbar® Magnetic Sink Strainer

Catalog No. F37787-0000 (blue) • F37787-0001 (green)



INSTRUCTIONS:

- Place strainer in sink drain in an upright position so that Scienceware® logo is visible.
- Flat edge allows strainer to be placed close to sink walls.
- Four strong magnets in each corner attract loose Spinbar® magnetic stirring bars and prevent them from falling into the drain.
- Periodically inspect and clean strainer as needed.



Help others choose the right product. Write a review at www.belart.com

661 Route 23 South, Wayne, NJ 07470

TEL: 973-694-0500 • FAX: 973-694-7199 • www.belart.com

937902113 – Rev. 2

Bel-Art – SP Scienceware assumes no obligation or liability for any advice furnished by it, or for results obtained with respect to these products. All such advice is given and accepted at the buyer's risk. © 2019. All Rights Reserved.

Bel-Art® – SP Scienceware® is a division of SP Industries, Inc.

