

Fluorescent Ruler

Catalog Number 13355-0001



The Fluorescent ruler can be in observing and photographing Ethidium Bromide stained nucleic acid agarose and acrylamide gels. The scale will fluoresce visibly upon illumination with ultraviolet light. This serves as a measure of nucleic acid mobility, and can be used to quantify the sizes of the nucleic acid fragments. This measurement will also allow the comparison of gels run at different times under the same conditions.



SP Bel-Art

2024 Broad Street, Pocomoke, MD 21851 USA

Tel: 973-694-0500 | cservice@scientificproducts.com

belart.com

Individual use varies. SP Bel-Art assumes no obligation or liability for any advice furnished by it or for results obtained with respect to these products. All advice is given and accepted at the user's risk. For warranty information and to register your product visit scientificproducts.com.

