

Aquet® Liquid Detergent

Catalog No. 17094-0030, 17094-0050

Dilute Aquet® to a 1% solution for general purposes. Although effective with cold water, a temperature of 90-100°F for Aquet® and rinse water can hasten drying.

Aquet® is a 25% solution in water of an aromatic polyglycol ether and is absolutely non-toxic and non-sensitizing to the skin. Safe, gentle and very effective wetting agent and detergent for laboratory glassware. It will not mar, etch or frost surfaces of lenses, cover glasses, cells, etc., and leaves no residue.



You May Also be Interested in:

Catalog No.	Description
17089-0000	Chromerge® chromic-sulfuric acid concentrate
17080-0200	Blowhard OS Extra® where total non-flammability is required



SP Bel-Art

2024 Broad Street, Pocomoke, MD 21851 USA

Tel: 973-694-0500 | cservice@scientificproducts.com

belart.com

Individual use varies. SP Bel-Art assumes no obligation or liability for any advice furnished by it or for results obtained with respect to these products. All advice is given and accepted at the user's risk. For warranty information and to register your product visit scientificproducts.com.

