


Chromatographic Column Holder

Catalog No. F18202-0000

Mount the holder on two support stands by sliding the support blocks onto the support stand rods and tightening the screws.

Loops move sideways in the base and can hold columns in between loops.

When using larger diameter columns, some loops may be removed to provide more space. This should be done during set-up to prevent dislodging of hanging columns.



SP Bel-Art

2024 Broad Street, Pocomoke, MD 21851 USA

Tel: 973-694-0500 | cservice@scientificproducts.com

belart.com

Individual use varies. SP Bel-Art assumes no obligation or liability for any advice furnished by it or for results obtained with respect to these products. All advice is given and accepted at the user's risk. For warranty information and to register your product visit scientificproducts.com.

