

Rotary Pipette Stand



Catalog No.18957-0000

Using the photograph as a guide, press the disc with concentric grooves onto one end of the short larger diameter tube, then press the disc with multiple holes onto the other end of the short, large diameter tube. Make sure these three parts are firmly pressed together

Slide this disc and tube assembly over the longer small diameter tube.

Press the open end of the longer tube with disc assembly into the base center.

All parts of this rack are made of polypropylene and autoclavable



SP Bel-Art

2024 Broad Street, Pocomoke, MD 21851 USA

Tel: 973-694-0500 | cservice@scientificproducts.com

belart.com

Individual use varies. SP Bel-Art assumes no obligation or liability for any advice furnished by it or for results obtained with respect to these products. All advice is given and accepted at the user's risk. For warranty information and to register your product visit scientificproducts.com.



918957001 – Rev 3