


Contact Plate Rack and Petri Dish Rack

Catalog No. F18990-0000, F18991-0000 & F18985-0100

Assembly:

1. Refer to drawing for orientation of parts.
2. Assemble the center divider to the base securing with one nut on each of the two threaded posts.
3. Assemble each of the rods to the base, the posts are keyed to fit with a specific orientation. Secure each rod with one nut on the threaded post.




You May Also be Interested in: Write-On™ Tape Dispensers Make labeling easy and fast

Producing legible labels is so much easier on a smooth, flat platform than on curved surfaces such as bottles or vials. The new Bel-Art Write-On™ Tape Dispenser provides this convenience by always having labeling tape positioned on the large writing platform, ready to make into a label. Select either the benchtop model with weighted base or economy model that can be wall mounted to save space. All models are easy to use by both left and right handed users.

Each dispenser can hold 5 or more tape rolls with widths from 12 - 25mm ($\frac{1}{2}$ - 1"). Ideal for use with the full line of Bel-Art, Write-On™ Tapes, F13463 series, available in 6 colors and three widths, 12, 19, and 25mm ($\frac{1}{2}$ ", $\frac{3}{4}$ " and 1"). This advanced tape adheres to any clean surface and resists moisture and most solvents. It withstands long exposure to temperatures as low as -73°C (-68°F) and up to 80°C (176°F) yet peels off without leaving a sticky residue.

Contact your preferred distributor to order these and other Bel-Art products.

Catalog No. F13472-0000 - Benchtop Model
F13473-0000 - Economy Model



Economy Model (shown)
Can be wall-mounted

Benchtop Model (not shown)
Includes weighted base



SP Bel-Art

2024 Broad Street, Pocomoke, MD 21851 USA

Tel: 973-694-0500 | cservice@scientificproducts.com

belart.com

Individual use varies. SP Bel-Art assumes no obligation or liability for any advice furnished by it or for results obtained with respect to these products. All advice is given and accepted at the user's risk. For warranty information and to register your product visit scientificproducts.com.



918990012 - Rev 6