


Velveteen Squares

Catalog No. H37848-0001 (Off White) & H37848-0002 (Black)

These 100% cotton, sterile Velveteen Replica Plating Squares are especially suited for use in conjunction with the SP Bel-Art Replica Plating Tool #F37848-0000.

The Velveteen Squares are supplied sterilized and ready for use. These may be re-used a number of times after careful washing and sterilization. To maintain sterility do not open the remaining bags until they are ready to be used. Seal tightly when not in use to reduce airborne contamination.




You May Also be Interested in:

The Bel-Blotter™ 96 Wells Replicating Tool

Replicate! Replicate! Replicate!

- Efficient, Cost Effective, Precise
- Capillary Action Ensures That Each Well Will Pick Up 10ml Per Tip
- 96 “Open Ended” Pipette Tips Mate Perfectly With Any 96 Well Plate
- Molded of Autoclaveable Polypropylene

Catalog No. F37876-0002



SP Bel-Art

2024 Broad Street, Pocomoke, MD 21851 USA

Tel: 973-694-0500 | cservice@scientificproducts.com

belart.com

Individual use varies. SP Bel-Art assumes no obligation or liability for any advice furnished by it or for results obtained with respect to these products. All advice is given and accepted at the user's risk. For warranty information and to register your product visit scientificproducts.com.

