


Mini LED Light Box

Catalog Number F37864-3000

- Remove battery cover
- Insert 6 AAA batteries (not included);
can also be powered by 120V AC adaptor
(not included)
- Use switch on lower left side to turn box on
- Gently place grid to cover entirety
of viewing area (if applicable)



You May Also be Interested in:



Extra Counting Grids (3)
F37864-3200



Colony Counter System
F37863-3000
complete system with
magnifier, colony counter,
and light box



Mini Magnifier
F37865-0000
for further enhanced
viewing



Colony Counter
F37862-0002
for keeping track of
colonies digitally



SP Bel-Art

2024 Broad Street, Pocomoke, MD 21851 USA

Tel: 973-694-0500 | cservice@scientificproducts.com

belart.com

Individual use varies. SP Bel-Art assumes no obligation or liability for any advice furnished by it or for results obtained with respect to these products. All advice is given and accepted at the user's risk. For warranty information and to register your product visit scientificproducts.com.



937864005 - Rev. 3