


Mouse Restrainer with Dorsal Access

Catalog No. F46401-0000

Instructions for Use:

This restrainer is designed to hold mice with body weights between 18 and 35 grams, and can be used for tail bleeds, tail tattooing, and for dorsal subcutaneous injections. Inspect restrainer for any signs of damage before using. Place the restrainer on a clean flat surface. If needed, the restrainer can be held in place with fasteners applied through the four slots located on the bottom of each corner. Do not leave the mouse unattended.



Your opinion counts! Write a product review at www.sp-belart.com



Bel-Art Products, a part of SP

661 Route 23 South, Wayne, NJ 07470 USA
Tel: 973-694-0500 | cservice.belart@spindustries.com
sp-belart.com

Individual use varies. Bel-Art Products/SP assumes no obligation or liability for any advice furnished by it or for results obtained with respect to these products. All advice is given and accepted at the user's risk. For warranty information visit belart.com/corporate.
For warranty information and to register your product visit sp-belart.com/corporate


946401006 – Rev 3

Mouse Restrainer with Dorsal Access

Catalog No. F46401-0000

Instructions for Use:

This restrainer is designed to hold mice with body weights between 18 and 35 grams, and can be used for tail bleeds, tail tattooing, and for dorsal subcutaneous injections. Inspect restrainer for any signs of damage before using. Place the restrainer on a clean flat surface. If needed, the restrainer can be held in place with fasteners applied through the four slots located on the bottom of each corner. Do not leave the mouse unattended.



Your opinion counts! Write a product review at www.sp-belart.com



Bel-Art Products, a part of SP

661 Route 23 South, Wayne, NJ 07470 USA
Tel: 973-694-0500 | cservice.belart@spindustries.com
sp-belart.com

Individual use varies. Bel-Art Products/SP assumes no obligation or liability for any advice furnished by it or for results obtained with respect to these products. All advice is given and accepted at the user's risk. For warranty information visit belart.com/corporate.
For warranty information and to register your product visit sp-belart.com/corporate

946401006 – Rev 3