

Super Strength, Dual Tested – Biohazard Disposal Bags

Meets ASTM Dart Impact and Tear Strength Requirements



- Meets 165g ASTM D1709 Dart Impact and 480g ASTM D1922 Tear Resistance tests required for Federal US DOT Interstate transport
- Meets all required regulations for use in Florida, Pennsylvania, and California
- Extra heavy-duty; 3 mil thickness
- Sterilization indicator patch darkens when steam sterilized for 20 minutes up to 135°C (275°F)
- Flat seam bags are made of high impact polypropylene

Catalog No.	Thickness	Dimensions (W x H)	Gallon Capacity	Quantity/Pack
F13167-1419	3.0 mil	36 x 48cm (14 x 19")	2 to 4	200
F13167-1923	3.0 mil	48 x 58cm (19 x 23")	5 to 9	200
F13167-2535	3.0 mil	64 x 89cm (25 x 35")	13 to 20	200
F13167-3138	3.0 mil	79 x 96cm (31 x 38")	20 to 30	200
F13167-3748	3.0 mil	94 x 122cm (37 x 48")	40 to 50	100