

Single Channel Timer

Three Key Operation for Ease of Use

- 99 minutes: 59 seconds Count Down/Count Up Function
- Last Time Set Memory
- Audible Alarm
- Stand and Magnetic Attachment



specifications

Timer Accuracy	± 5 seconds/day
Display Resolution	1 second
Display Size	48.3W x 15.2mmH (1.9 x 0.6")
Battery	One 1.5 volt, type AAA or equivalent (included)
Product Size	61W x 61H x 15.2mmD (2.4 x 2.4 x 0.6")

BEL-ART CATALOG NO.	DESCRIPTION	CASE QTY
B61700-3400	Single Channel Timer	10

Ten Button Timer

Ten Buttons Allow Direct Input

- Count Down/Count Up Function
- Audible and Light Alarms
- Pocket Clip, Foldaway Stand and Magnetic Attachment



specifications

Timer Accuracy	± 5 seconds/day
Display Resolution	1 second
Display Size	53.3W x 22.9mmH (2.1 x 0.9")
Battery	Two 1.5 volt, type AAA or equivalent (included)
Product Size	81.3W x 71.1H x 30.5mmD (3.2 x 2.8 x 1.2")

BEL-ART CATALOG NO.	DESCRIPTION	CASE QTY
B61700-3500	Ten Button Timer	10

Four Channel Timer with White Board

Track Time With Built-In White Board

- Four Channels
- Counts Up To or Down From 23 hours: 59 minutes and Can Be Set In Increments of 1 Minute
- Audible Alarms
- Key Press Confirmation
- Foldaway Stand or Magnetic Attachment



specifications

Timer Accuracy	± 5 seconds/day
Display Resolution	1 second
Display Size	43.2W x 53.3mmH (1.7 x 2.1")
White Board Size	61W x 68.9mmH (2.4 x 2.7")
Battery	Two 1.5 volt, type AAA or equivalent (included)
Product Size	114.3W x 88.9H x 25.4mmD (4.9 x 3.5 x 1")

BEL-ART CATALOG NO.	DESCRIPTION	CASE QTY
B61700-3700	Four Channel Timer with White Board	10

All digital timers and stopwatches feature a **Certificate of Calibration** and individual serial number, and are tested and calibrated in H-B's exclusive triple accredited/registered ISO/IEC 17025:2017, A2LA accredited laboratory (Cert #2448.01), against equipment whose calibration is traceable to NIST. CE marked, RoHS compliant.



www.belart.com • cservice@scientificproducts.com • +1-800.220.5171

© 2025 SP Bel-Art. All rights reserved.
SP (scientificproducts.com) is a manufacturer of biopharma processing equipment, life sciences equipment, scientific glassware and laboratory supplies serving pharmaceutical, biotechnology, educational, food and beverage and industrial markets.
Care has been taken to ensure that the information in this literature is accurate, however we do not represent or warrant that the information is current or accurate. In some cases, the product dimensions may be rounded for presentation. If specifications are critical to your application please contact cservice@scientificproducts.com to confirm the data from this catalog.
All trademarks referenced in this document are registered or common law trademarks of Bel-Art Products, Inc. or its affiliates in the U.S. and/or other countries.



PRECISION INSTRUMENTS FOR
Temperature | Density | Time | Calibration



SP-01006-RevC

Triple Alarm Single Channel Timer

Perfect for High Noise Locations

- Vibration, Light and Loud Buzzer Sound Alarms
- 99 minutes 59 seconds count up/down
- Key Press Confirmation
- Pocket Clip, Foldaway Stand and Magnetic Attachment
- Clock function included



specifications

Timer Accuracy	± 5 seconds per day
Display Resolution	1 second
Display Size	35 x 20mmH (1.4 x 0.8")
Battery	Two 1.5 volt, type AAA or equivalent (included)
Product Size	56W x 75H x 32mmD (2.2 x 3 x 1.3") exclude clip portion

BEL-ART CATALOG NO.	DESCRIPTION	CASE QTY
B61700-3600	Triple Alarm Single Channel Timer	10



B61700-5300



B61700-5200



B61700-5000



B61700-5100

Timers



B61700-2500



B61700-2600



B61700-3000



B61700-3200



B61700-3300



B61700-2900

BEL-ART CATALOG NO.	CHANNELS	HOURS	MINUTES	SECONDS	MEMORY	CLOCK	BATTERY	MOUNT	DIMENSIONS	ACCURACY	CASE QTY
DURAC® Timers											
B61700-2500	1	-	60	-	NO	NO	-	-	57W x 57mmH (2.2" x 2.2")	± 3 min	20pcs
B61700-2600	1	-	99	59	YES	NO	(1) AAA	Fold away stand, pocket clip, magnetic back	67W x 70mmH (2.6" x 2.8")	± 3 seconds	20pcs
B61700-3000	1	23	59	59	YES	YES	(1) LR44	Fold away stand, pocket clip, magnetic back	70W x 51mmH (2.8" x 2")	± 3 seconds	10pcs
B61700-3200	3	19	59	59	YES	YES	(2) LR44	Fold away stand, pocket clip, magnetic back	62W x 118mmH (2.5" x 4.6")	± 3 seconds	20pcs
B61700-3300	4	99	59	59	YES	YES	(1) AAA	Fold away stand, pocket clip, magnetic back	80W x 80mmH (3.1" x 3.1")	± 3 seconds	20pcs
B61700-2900	1	19	99	59	YES	NO	(1) AAA	Fold away stand, pocket clip, magnetic back	64W x 64mmH (2.5" x 2.5")	± 5 seconds	20pcs

Stopwatches

BEL-ART CATALOG NO.	TYPE	TIME	CHRONOGRAPH	SPLIT/LAP TIME	PACER	MEMORY	CLOCK	ALARM	CALENDAR	CONSTRUCTION	BATTERY	DIMENSIONS	ACCURACY	CASE QTY
DURAC® Stopwatches														
B61700-5300	Digital	-	1/100 SEC	-	-	1	YES	YES	YES	Plastic	(1) CR2032	63W x 80mmH (2.5" x 3.1")	0.01 sec	20pcs
B61700-5200	Digital	10HR	-	99	10/320 SEC	30	YES	NO	YES	Plastic	(1) CR2032	63W x 83mmH (2.5" x 3.3")	0.01 sec	10pcs
B61700-5000	Analog	15MIN 1/10SEC Intervals	-	-	-	NO	NO	NO	NO	Copper Chromium Plated	-	63W x 83mmH (2.5" x 3.3")	0.2 sec	5pcs
B61700-5100	Analog	30MIN 1/5SEC Intervals	-	-	-	NO	NO	NO	NO	Copper Chromium Plated	-	63W x 83mmH (2.5" x 3.3")	0.2 sec	5pcs

All digital timers and stopwatches feature a **Certificate of Calibration** and individual serial number, and are tested and calibrated in H-B's exclusive triple accredited/registered ISO/IEC 17025:2017, A2LA accredited laboratory (Cert #2448.01), against equipment whose calibration is traceable to NIST. CE marked, RoHS compliant.