

Secador® 5.0 Desiccator Cabinets

The Secador 5.0 is available in three models. A non-electric model for use with desiccant, two desiccant free electric models, and a digitally adjustable model for easy control of relative humidity. These desiccators are an ideal solution for storage of laboratory reagents and samples, optical equipment, electronics, and 3D printing filament.



- High density foam gaskets keep ambient humidity out and dry air inside
- Clear polycarbonate panels are durable and allow for easy viewing of contents
- LCD displays internal temperature and humidity level
- Non-electronic model can be used with any Bel-Art desiccant product
- Includes 3 adjustable and removable shelves

Catalog No.	Description	Electronic Dehumidification	Adjustable Humidity Control	Humidity Control Range	Internal Dimensions (W x D x H)	External Dimensions (W x D x H)
H42057-0001	Secador 5.0 Desiccator Cabinet	No	No	Desiccant required	28cm x 34cm x 50cm	30cm x 38cm x 52cm
H42057-1003	Secador 5.0 Auto-Desiccator Cabinet	Yes 110-240VAC	No	Maintains 25 to 40%	28cm x 34cm x 50cm	30cm x 38cm x 63.5cm
F42057-1004	Secador 5.0 Adjustable Auto-Desiccator Cabinet	Yes 110-240VAC	Yes	Adjustable from 20 to 40% rH	28cm x 34cm x 50cm	30cm x 38cm x 63.5cm
F42057-1001	Secador 5.0 Shelf					

Desiccant products for use with H42057-0001 Secador 5.0 Desiccator Cabinet

- F42045-0350 – Desiccant Bag, 350g
- F42045-0500 – Desiccant Bag, 500g
- F42045-1000 – Desiccant Bag, 1000g
- F42046-0000 – Desi Can Reusable Desiccant Canister
- F42049-0100 – Reusable Silica Gel Desiccant Cartridges, 72mm