



BEL-ART PRODUCTS

Press Release
October 2, 2014

For Immediate Release
Contact:

Teri L. Rocco
trocco@belart.com – 1-800-4-BELART
Bel-Art Products, Inc.

Squeeze Rite™ Wash Bottle Helps Reduce Waste in the Lab

Wayne, NJ – October 2, 2014 – At a time when many are feeling their laboratory budgets being squeezed, Bel-Art Products has introduced a new kind of squeeze aimed at containing lab costs by reducing waste. New to market is the Scienceware® Squeeze Rite™ Wash Bottle, catalog no. F11717-0030. Innovative in shape and design, the dip tube position within the bottle reduces the waste of reagents, solutions and other liquids by allowing a more complete dispensing of those liquids.



“Easier to squeeze when compared to standard 500ml bottles, its unique shape fits the hand comfortably, and can help reduce hand fatigue,” according to Edmund Murphy, Category Manager for Bel-Art Products.

Made of low-density polyethylene (LDPE), the Squeeze Rite™ Wash Bottle has a generous, 53mm diameter mouth that makes it a cinch to fill from larger bottles and other primary containers. Labeling is simple using a marker or label compatible with LDPE.

Suggested list price for the Scienceware® Squeeze Rite™ Wash Bottle, catalog no. F11717-0030, is \$52.00 for a pack of 6.

For more information on the Scienceware® Squeeze Rite™ Wash Bottle, or the complete line of laboratory offerings from Bel-Art Products, call 1-800-4BELART, email info@belart.com or visit www.belart.com. To request a FREE copy of the Bel-Art Products 2011 Catalog, [click here](#).

661 Route 23 South, Wayne NJ 07470
Tel: 973-694-0500 • Fax: 973-694-7199 • email: info@belart.com • www.belart.com

Subsidiaries: Maddak, Inc. • Nutech Manufacturing Corporation • Applied Coatings, Inc.



BEL-ART PRODUCTS

Bel-Art Products simplifies tasks with innovative product solutions for science, industry and healthcare. Over 4500 items from safety wash bottles to fume hoods are manufactured and marketed under well-known brand names including Scienceware® laboratory products, Spinbar® magnetic stirring bars, Magic Touch™ icewares, Sterileware® sampling tools, H-B Instrument™ thermometers and hydrometers, Secador® desiccator cabinets, Flowmi™ Cell Strainers, and Poxymat® racks. Established in 1946, Bel-Art Products is an SP Industries company headquartered in Wayne, NJ.

#

661 Route 23 South, Wayne NJ 07470
Tel: 973-694-0500 • Fax: 973-694-7199 • email: info@belart.com • www.belart.com

Subsidiaries: Maddak, Inc. • Nutech Manufacturing Corporation • Applied Coatings, Inc.