

## Lab-Aire® II Drying Racks

Laboratory space is always in high demand making the small footprint of the Lab-Aire® II Drying Racks a practical solution for drying and storing fragile, frequently used labware.

Choose from **Non-electric Benchtop** or **Wallmount** drying racks or accelerate drying with the warm, flowing air from **Electric Benchtop** models. All styles have removable pegs in 3 different lengths that can be positioned for optimal use. Extra accessories and a variety of different sizes/configurations let you create a drying system that works for you.

### Wallmount - Non-electric

- Space-saving, move drying off the bench
- Easily mounts on walls or standard pegboard
- Single-sided; can be custom made with additional height

### Benchtop - Non-electric

- Portable, use anywhere convenience
- Offers the most choices in size
- Single- and Double-sided models

### Benchtop - Electric

- Accelerate drying with warm, circulating air
- Use with or without power cycle
- Single- and Double-sided models: 120V or 230V



**Wallmount - Non-electric**  
Single-sided, 2 tier



**Benchtop - Non-electric**  
Double-sided, 1 tier



**Benchtop - Electric**  
Single-sided, 3 tier

### Increase Drying Space Three Ways

1. Go up in height – 1, 2 and 3 tier models available in all styles
2. Double-sided benchtop models double your capacity with minimal increase in footprint
3. Use multiple Lab Aire® II Dryers

### Increase Flexibility

Each Lab Aire® II Dryer comes with an assortment of pegs in various lengths. Electric models also have baskets for holding small labware. Additional pegs and baskets\* to customize your dryer can be purchased separately.

\* Baskets fit electric models only.

## Lab-Aire® II Dryers Quick Comparison Chart

Select from 12 different models with up to 114 peg positions.

- Rugged ABS plastic will never corrode and easily wipes clean
- Drip/storage trays with scalloped wells hold long thin labware
- Small (3"), Medium (4"), and Large (5") pegs can be moved and repositioned to hold a variety of labware including beakers, flasks, bottles, test tubes and cylinders
- Extra pegs and drying baskets available

	Style	Wallmount			Dimensions (W x D x H)	Peg Slots	Peg Sizes*	Model
			Benchtop	Electric				
Wallmount Non-Electric	Single 1 Tier	X	X		37.5 x 12.7 x 24cm (14.75 x 5 x 9.4")	19	S - 6, M - 7, L - 6	18933-0011
	Single 2 Tiers	X			37.5 x 12.7 x 41.7cm (14.75 x 5 x 16.4")	38	S - 12, M - 14, L - 12	18933-0014
	Single 3 Tiers	X			37.5 x 12.7 x 59.4cm (14.75 x 5 x 23.4")	57	S - 18, M - 21, L - 18	18933-0013
Benchtop Non-Electric	Single 1 Tier	X	X		37.5 x 12.7 x 24cm (14.75 x 5 x 9.4")	19	S - 6, M - 7, L - 6	18933-0011
	Single 2 Tiers		X		37.5 x 17.8 x 39cm (14.75 x 7 x 15.4")	38	S - 12, M - 14, L - 12	18933-0012
	Single 3 Tiers		X		37.5 x 17.8 x 57cm (14.75 x 7 x 22.4")	57	S - 18, M - 21, L - 18	18933-0015
	Double 1 Tier		X		37.5 x 25.4 x 21.3cm (14.75 x 10 x 8.4")	38	S - 12, M - 14, L - 12	18933-0021
	Double 2 Tiers		X		37.5 x 25.4 x 39cm (14.75 x 10 x 15.4")	76	S - 24, M - 28, L - 24	18933-0022
	Double 3 Tiers		X		37.5 x 25.4 x 57cm (14.75 x 10 x 22.4")	114	S - 36, M - 42, L - 36	18933-0023
Benchtop Electric	Single 2 Tiers	X		X	42.5 x 19 x 40cm (16.75 x 7.5 x 15.7")	19	S - 6, M - 7 L - 6 2 Baskets, 7 Peg Plugs	18819-0012 (120V) 18819-2301 (230V)
	Single 3 Tiers	X		X	42.5 x 19 x 57.7cm (16.75 x 7.5 x 22.7")	38	S - 12, M - 14 L - 12 2 Baskets, 14 Peg Plugs	18819-0013 (120V) 18819-2302 (230V)
	Double 2 Tiers	X		X	42.5 x 25.4 x 40cm (16.75 x 10 x 15.7")	38	S - 12, M - 14 L - 12 4 Baskets, 14 Peg Plugs	18819-0022 (120V) 18819-2303 (230V)
	Double 3 Tiers	X		X	42.5 x 25.4 x 57.7cm (16.75 x 10 x 22.7")	76	S - 24, M - 28 L - 24 4 Baskets, 28 Peg Plugs	18819-0023 (120V) 18819-2304 (230V)

\*S - Small • M - Medium • L - Large